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La FDA autorise GentleWaves **29**

FDA Clears GentleWaves

American Society for Laser Medicine and Surgery **33**

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Mediform Group extends its activity in Europe

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Read precise tips on how to inject the endoprosthesis Bio Alcamid™ safely by using the "Volume Technique" but also to do additional injections whenever it is needed.

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Learn about the original "Implantation in X Technique" with long and barbed suspension threads, for dramatic aesthetic results.

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Radiofrequency device with a unipolar and bipolar double configuration, designed for face, body and neck remodelling as well as the treatment of cellulite.

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First ablative and non ablative system combining laser, UVB and IPL technologies on a single platform for a wide spectrum of clinical indications.

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VIII INTERNATIONAL SYMPOSIUM
On Laser, Intense Pulsed Light and other Dermocosmetic Techniques

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Neck, oval, cheeks (in association with Exopeel) after 4 months.

Aptos threads Implantation in X

Doctor Denis Guillo,
Aesthetic Doctor in Lyon (France)

Foreword

The APTOS™ threads need no presentation – barbed threads, known as Russian threads, – which since 1997, and thanks to their inventor, Doctor Marian Sulamanidze, have shaken the small world of non-surgical lifting.

The new reflection resulting from this technique is interesting, and seems a little puzzling to certain doctors.

Indeed, we are witnessing a new and unique competitor of the surgical lift that integrates the notions of the recovery potential of stretched skin, as well as the time required and that it is conservative; which makes it difficult to injure this very valuable tissue.

Even if the basic idea is very attractive, we think that the results obtained with short threads may be improved.


The evolution: Appearance of the X.

The improvement applies to two "details" that make it possible to make up for the shortcomings of the initial short threads.

1. The length, precisely. It must be increased for the sagging area to be linked directly to the fixation area. Threads of 23 cm are available, but that is still insufficient due to the maximum needle length (15 cm).
2. Securing of the paired threads by a central knot makes the implantation unit an X. This recent evolution, however small it may seem, has major advantages.
 - The traction permitted on the threads is very powerful,
 - The implantation length of the threads is maximal, particularly in the scalp where it was really necessary, which allowed the thread to migrate towards the bottom even although they were well balanced.
 - Migration is now no longer possible because the fibrosis blocks the knot.
 - The reversibility of the operation, which is very reassuring for patients, is rendered very easy by the location of the knot.


Implantation technique on neck, jowls and cheekbones.

The eyebrows, a lot more difficult, are not affected by this type of implantation. They will be addressed in a forthcoming presentation of a new V implantation, which will make it possible to forego about the temporary lifting and its scars.

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
Schema marquage.
Preoperative markings.

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
Implanting an X is not as easy as implanting a single thread, and the down side of these advantages is the greater technical complexity.


1. *The marking, very important, requires a certain degree of deliberation and makes it possible to obtain a personalised implantation outline. It is performed with the Onyx red marker and cleaned with acetone, with the patient seated – cf. diagram–.*
2. *Anaesthesia is with xylocaine with 2% adrenaline– 20% bicarbonate all along the thread line. This is the unpleasant part.*
3. *Implantation. The entry point is zygomatic and pre-mastoid. It is perforated with a Nokor needle. The 18 G double bevel needle, bent, handled with a needle holder, goes from the centre and moves towards the base of the ptosis. It comes out again after the corner of the mouth and the nasolabial*

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
fold. Insert one of the bottom branches of the X into the temporary tunnel formed by the needle, then extract the needle through the external end. Do the same for the other branch, then the two high branches that are implanted on the scalp. The knot will be inserted automatically.

4. **Fixation:** finally, with the patient seated, the threads are pulled upwards, balanced and then cut immediately.

The combination of the threads with other techniques – volume restructuring by lipofilling or smoothing by peeling– is perfectly possible extemporarily or secondarily, as well as a secondary surgical lifting if the result is unsatisfactory.

Follow-UP

Pains and ecchymoses –rare haematomas– predominate and may last 10 to 15 days. In the early days, a sometimes important triangulation of the face with refinement of the chin and elongation of the temples –due to the raised tissues, anaesthesia, oedemas– disturbs patients quite a lot.

You may suggest that patients do not socialise for a few days. Generally speaking, results are good, with a very good patient return.

The folds are visible –since there is no skin resection– in the preauricular area for one to two weeks and are concealed by the hair.

All independent observers –particularly hairdressers– find that results are even better after 3 - 4 months.

Conclusion

This variant of the barbed thread implantation technique, powerful and subtle, seems to be bringing long barbed threads up to the rank of an efficacious non-surgical lifting solution. Some time may need to elapse before this technique is fully appreciated, but never before have such impressive and immediate results been obtained with short threads. It is obvious that further progress, particularly with regard to thread quality, may be expected.

This treatment has far fewer limitations than the surgical act, which renders it feasible for a very large number of patients, and while this solution cannot yet replace the lifting in cases of major sagging, it may be reserved for incipient sagging with bad indication for surgical lifting.

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